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
Microscope Digital Camera M BASE SERIES

EN	User Manual
CZ	Návod k použití
DE	Bedienungsanleitung
ES	Guía del usuario
PL	Instrukcja obsługi
RU	Инструкция по эксплуатации

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Microscope Eyepiece Digital Cameras Levenhuk M BASE Series

Never look at the sources of bright light or lasers through the camera or the microscope: it will cause PERMANENT DAMAGE TO YOUR EYES!
--

Levenhuk eyepiece digital cameras have been designed exclusively for microscopes. They can be installed on any type of optical microscopes, be it biological or metallurgical, mono or stereo microscope. The observed view may be transferred directly to your PC, and with the software included in the kit you can view and edit the resulting image.

Features

- Levenhuk M BASE Series includes cameras with resolutions ranging from 0.3 to 5 megapixels;
- The camera allows you to watch images on the screen in a real-time mode, and to capture photos and videos;
- The camera renders realistic colors with high level of detail;
- Levenhuk M BASE Series cameras are compatible with microscopes of 23.2mm eyepiece tube diameter (the camera is installed in place of the eyepiece);
- The standard kit includes Levenhuk software necessary to view and edit your images (brightness, image size, etc.);
- The camera can be charged and connected to your PC with USB-cable.

Package

Levenhuk digital camera, USB cable, installation CD (setup wizard, drivers, software), user manual and warranty.

System requirements

Mac OS 10.12, Linux Ubuntu 14.04, Windows Vista/7/8/10 (32-bit or 64-bit), CPU IntelCore 2 or above; RAM 2GB, USB port 2.0, CD-ROM.

Caution

- Never look at the sources of bright light or lasers through the camera or the microscope: it will cause PERMANENT DAMAGE TO YOUR EYES!**
- Do not disassemble the camera on your own.
- Keep the camera away from moisture; do not use it in rainy weather.
- Keep the camera away from shock or excessive pressure. Do not overtighten the locking screws.
- Keep the camera away from hazardous environment, home and car heaters, incandescent lamps or open fire.
- When cleaning the lenses, first blow any dust or debris off the surface or wipe them off with a soft brush. Then wipe the lens with a soft cleaning cloth slightly damp with spirits or ether.
- Seek medical advice immediately if a small part or a battery is swallowed.**

Levenhuk International Lifetime Warranty

All Levenhuk telescopes, microscopes, binoculars and other optical products, except for accessories, carry a **lifetime warranty** against defects in materials and workmanship. Lifetime warranty is a guarantee on the lifetime of the product on the market. All Levenhuk accessories are warranted to be free of defects in materials and workmanship for **six months** from date of retail purchase. The warranty entitles you to free repair or replacement of the Levenhuk product in any country where a Levenhuk office is located if all warranty conditions are met. For further details please visit our web site: https://www.levenhuk.com/warranty. If warranty problems arise, or if you need assistance in using your product, contact the local Levenhuk branch.

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Levenhuk® is a registered trademark of Levenhuk, Inc.
© 2006-2017 Levenhuk, Inc. All rights reserved.
Levenhuk, Inc. 924-D East 124th Ave. Tampa, FL 33612 USA